

H2D2

To calibrate

With SynerJY program open

Collect

Experiment setup

Click on "Detectors" icon

Start: 546,1 nm

End: 546.3 nm

Increment: 0.001 nm

Integration time: 0.01 seconds

High voltage supply: input 750 volts

Click on/off box to green

Click on "Monos" icon

Slits

Front entrance: 0.02nm

Front exit: 0.02 nm

Turn on Hg source

Position at entrance

Manual slits set at 1

Click Run in SynerJY program

When run completes

Determine the center wavelength between the two observed peaks

The listed wavelength in a vacuum for this peak of mercury is 5460.735 Angstroms

from <http://physics.nist.gov/RhysRefData/Handbook?Tables/mercurytable2.htm>

Log on to: <http://emtolbox.nist.gov/wavelength/edlen.asp>

Input for vacuum wavelength: 546.0735 nm

Input current air temperature

Input current atmospheric pressure

Input current air humidity

Hit: Calculate Wavelength in Ambient Air and Refractive Index of Air

In SynerJY

Move the grating to the center position

Select to "Monos" button

Click the "Calibrate" button

In the Calibration Dialog window

Current position should read the value for the determined center point

In the Calibrate Position box, input the determined center point wavelength

Click OK

The monochromator is now calibrated for the wavelength of Mercury in Air at current conditions.